

	Appearance	Assay	Assay	Assay	Assay	Moisture (KF)	Hardness	Dissolution	Dissolution	Dissolution
	Appearance	Aspirin	Impurity A	Impurity B	Total Impurities	Moisture (KF)	Hardness	Aspirin, individual @ 45 min	Aspirin, mean @ 45 min	Aspirin, rsd @ 45 min
		% label claim	%	%	%	%	kp	% I.c.	% I.c.	%
1000										
Initial	White to off white tablet	100.0	0.1	0.1	0.2	1.1	90.0	99, 98, 97, 99, 98, 97	98	0.1
25°C ± 2°C/60% RH ± 5% RH										
1 month	White to off white tablet	98.9	0.1	0.2	0.3	1.2	90.5	99, 98, 97, 98, 97, 99	98	0.1
3 month	White to off white tablet	99.1	0.2	0.1	0.3	1.2	90.7	98, 99, 97, 99, 98, 97	98	0.1
6 month	White to off white tablet	98.4	0.2	0.2	0.4	1.3	90.1	98, 99, 99, 99, 99, 98	99	0.1
9 month	White to off white tablet	97.6	0.3	0.3	0.6	1.5	89.1	97, 97, 98, 97, 98, 97	98	0.1
12 month	White to off white tablet	96.5	1.0	0.7	1.7	1.8	87.1	97 (U), 98 (U), 96 (U), 98 (U), 97 (U), 96 (U)	98 (U)	0.1 (U)
40°C ± 2°C/75% RH ± 5% RH										
12 month	-----	-----	-----	-----	-----	-----	-----	79 (U)	-----	-----